


# AIRSTRIP®

POCKET GUIDE

**AIRSTRIP ONE®**  
for obstetrics



# INTRODUCTION

**AirStrip ONE for Obstetrics delivers visibility into monitored patient data in the following categories:**

- Real-time monitor viewing
- Patient Demographics

**Supported data include:**

- Waveforms
- Dilatation Status and Cervical Exams
- Annotations
- Maternal Vitals

**Launching and using AirStrip ONE:**

1. Tap the AirStrip ONE app
2. Tap the site's name
3. Log in with your site credentials


**To get started or get  
help, call 877-258-5869  
[www.airstrip.com](http://www.airstrip.com)**

**Important Note: Users must  
contact AirStrip at (877) 258-5869  
for re-installation of AirStrip  
ONE if they acquire a new mobile  
device.**

**For complete instructions please  
see the user guide.**

# ACCESSING PATIENT DATA

Locate a patient using SEARCH, or from the default patient list.



# PATIENT CENSUS SCREEN

A screenshot of the Patient Census screen. The screen shows a list of patients with various details. Red dashed lines and boxes connect text annotations to specific elements on the screen. The annotations include: 'Tap the Navigation Menu button to view the Navigation List and to Search for a patient.' pointing to the navigation menu icon; 'Unit name' pointing to the 'North Wing' header; 'Dilation Status' pointing to a green dot; 'Estimated Gestational Age(EGA)/Effacement/Station and Time of last exam' pointing to a date and time; 'Membrane Status' pointing to a green dot; 'To sort or filter the patient list, select the Sort or Filter button.' pointing to the sort and filter icons; 'Doctor Last Name Nurse Last Name Patient Location' pointing to the patient details; 'Note: Display of the doctor and nurse's name is configurable by Site.' pointing to the patient details; 'Tap to view the Patient Summary screen.' pointing to the right arrow; and 'Scroll to view additional patients.' pointing to the bottom of the list.

Tap the Navigation Menu button to view the Navigation List and to Search for a patient.

Unit name

Dilation Status

Estimated Gestational Age(EGA)/Effacement/Station and Time of last exam

Membrane Status

To sort or filter the patient list, select the Sort or Filter button.

Doctor Last Name  
Nurse Last Name  
Patient Location  
**Note:** Display of the doctor and nurse's name is configurable by Site.

Tap to view the Patient Summary screen.

Scroll to view additional patients.

# DILATION STATUS

Patients that are 1-4.9 centimeters dilated



Ashford, Allison C.

Intact

37.1 • G2P2 • 3/40/1 @ 01:04

Obst (MC)

Long

LDR-4

Patients that are 5-7.9 centimeters dilated



Collinsworth, Gena

Ruptured

39.1 • G4P2 • 6/90/1 @ 13:20

Obst (MC)

Long

LDR-7

Patients not yet documented



Helms, Natalie

Intact | Epidural

36.1 • G1P0 • -- / -- / --

Obst (MC)

Long

LDR-4

Patients that are 8-10 centimeters dilated



Wallace, Georgina

Ruptured | Epidural

39.2 • G3P2 • 8/100/2 @ 23:01

Obst (MC)

Long

LDR-9

# PATIENT SUMMARY SCREEN

The header displays the patient name, age, gender, hospital name, unit, and bed number.

Tap the **Action Overflow** button to share the screen with another AirStrip ONE user, toggle strip colors or change column width.

Dilation and Membrane status.

Tap the X to close the screen.

EGA/G, P/Dilatation/Effacement/Station.

Tap the caret button to show/hide the waveforms.

Pinch/zoom to zoom in/out on strips. Note the change in time and % of the strip chart displayed.

Tap the **Live (real-time)** button to return to live (real-time) view.

Swipe right view history for up to 24 hours.

Tap the **Annotation** icon to view an annotation.

Drag the screen divider up to view tiles.

**Time Ago** indicator displays when accessing historical data.

Tap a tile to view more patient data.



# MATERNAL VITAL SIGNS SUMMARY TILE



High - Low range and most recent Vitals for up to 23 hours for:

- T (Temperature)
- HR (Heart Rate)
- SBP (Systolic Blood Pressure)
- DBP (Diastolic Blood Pressure)
- RR (Respiration Rate)
- O2% (Oxygen Saturation)

Note: Summary tiles display a summary of the patient data available on the expanded tile. Once a summary tile is tapped, the expanded tile displays.

# AIRSTRIP®

[www.airstrip.com](http://www.airstrip.com)

# AIRSTRIP®

2915 W. Bitters Rd., Suite 215  
San Antonio, TX 78248  
+1 210 805 0444  
www.airstrip.com

EU Authorized Representative;  
Emergo Europe B.V.  
Prinsessegracht 20  
2514 AP The Hague, The Netherlands  
+31 70 345 8570



Cardioscan Services Pty Ltd.  
293 Camberwell Rd., Suite 301  
Victoria, 3124 Australia  
+61 3 9832 2222  
0 3 9832-2222

© 2017-2018. AirStrip Technologies. All rights reserved. iPhone® and iPad® are registered trademarks of Apple Inc. Trademarks not belonging to AirStrip Technologies are the property of their respective companies. Screen images contained within are for illustrative purposes. Actual screens may vary by implementation, site and user.  
AirStrip ONE for obstetrics iOS pocket guide V 1.7 2018-10-11



AirStrip Technologies, Inc.